

PRODUCT CODE: CS3039

PRODUCT NAME: CleanSpace™ Steri-Plus Exhalation Valve Filter
(pk 20)



Description

The CleanSpace Steri-Plus Exhalation Filter is suitable for use with the CleanSpace HALO. The 'Steri-Plus' (CS3038 and CS3039) does not contribute to the system's respiratory protection and must be used alongside either a CleanSpace HALO Particulate Filter (CS3022) or a CleanSpace HALO Bio Particulate Filter (CS3025). The Steri-Plus is designed to filter the air exhaled by the wearer and consequently protects those not wearing a respiratory protection system. The Steri-Plus Exhalation Filter should be changed if it has been contaminated with organic solvents (for instance isopropanol, xylene or toluene), if it has become visibly soiled or if it has passed its expiry date.

Note: The Steri-Plus shall only be used in combination with a HALO Half Mask fitted with the PAF-1111 Exhalation Valve Cover. The Steri-Plus shall NOT be used with the CS3024 BioHood.

Test Specification

The Steri-Plus Exhalation Filter is batch tested to ensure its effectiveness. This batch testing ensures that at a flow rate of 30 LPM the penetration through the filter is less than 1% and the exhalation resistance is less than 15 mm H₂O.

Features

- Used with the revolutionary CleanSpace HALO - A light weight PAPR with no hoses or belts.
- Materials: Nonwoven electrostatic media.
- Quickly and easily fitted and removed from the Exhalation Valve Filter Case (CS3038).
- Sold in a pack of 20 (CS3039).

Specifications and materials

- Weight: 750g (pkt 20) Package Dimensions: 200mm x 270mm
- Materials: Nonwoven electrostatic media, polypropylene.
- Storage and Use: -10°C to +55°C (-4°F to +131°F) at <90% relative humidity. Store away from direct sunlight, grease and oil.
- Cleaning: Do not clean, dispose of when a change of Exhalation Filter is required.
- Only to be used alongside a CleanSpace Particulate Filter. The Steri-Plus does not contribute to the system's respiratory protection.

Suitable Applications

Healthcare, pharmaceutical production, research, diagnostic laboratories and emergency responders. Suitable for protection against particulates including airborne biohazards.

Training

Online training available with verification for compliance purposes.
Contact sales@cleanspacetechnology.com

Limitations

CleanSpace respirators are air filtering, fan assisted positive pressure masks and designed to be worn in environments where there is sufficient oxygen to breathe safely. Do not use the CleanSpace in IDLH atmospheres, to protect against gases/vapours that cannot be filtered, or in Oxygen enriched or deficient atmospheres.