



FIT TESTING - CLEANSPACE HALF MASK

When conducting a quantitative fit test, CleanSpace Half Masks are required to be assembled with the following items:

- CleanSpace Half Mask Fit Testing Adaptor (PAF-0025).
- Two (2) CleanSpace Grommets (PAF-1036 for a pack of 6).

CST1043 – CleanSpace Half Mask Fit Test Kit includes 1 x CleanSpace Half Mask Fit Testing Adaptor, 6 x CleanSpace Grommets + Instructions

CleanSpace fit test grommets are inserted into both arms of the Half Mask. The black circular diaphragm of the grommet should be oriented towards the face. The blades should be fully inserted into the plastic oval openings at the end of the Half Mask arms.

NOTE: Fit test grommets MUST be removed when a CleanSpace respirator is in operational use (PAPR- Power ON/ motor supplying positive airflow throughout the respirator)



PREPARE THE CLEANSPACE HALF MASK FOR FIT TESTING		
<p>1</p>	<p>2</p>	<p>3</p>
<p>Clean the Half Mask before each use.</p>	<p>Insert a grommet into each arm of the Half Mask.</p>	<p>Remove exhalation valve assembly (with cap).</p>
<p>4</p>	<p>5</p>	<p>6</p>
<p>Attach Half Mask adaptor to exhalation valve.</p>	<p>Connect clear fit testing tube.</p>	<p>Re-attach Half Mask adaptor into the mask.</p>

Follow the procedure as outlined in the Manuals provided with the Fit Testing Equipment.

