



## CleanSpace™ HALO Cleaning and Storage Plug



<b>PRODUCT CODE</b>	CS3011
<b>PRODUCT NAME</b>	CleanSpace Cleaning and Storage Plug
<b>DESCRIPTION</b>	<p>The CleanSpace HALO Cleaning and Storage Plug is designed to be used with the CleanSpace HALO (CS3001) power unit and the CleanSpace WORK (CS3001WORK) power unit.</p> <p>This accessory is used to seal both bellows and filter air inlet. Please remove the filter from the CleanSpace power unit before use.</p> <p>During cleaning, all openings must be sealed. It is also recommended that all openings are sealed after cleaning and decontamination and during storage.</p>
<b>FEATURES</b>	<ul style="list-style-type: none"><li>• Clips easily into position</li><li>• Facilitates a high level of cleaning</li><li>• Keeps dust and moisture from entering during and after cleaning, decontamination, and storage</li></ul>
<b>SPECIFICATIONS</b>	<ul style="list-style-type: none"><li>• Product weight: 3 oz (85 g)</li><li>• Product dimensions: 170 mm x 50 mm x 40 mm (approximately)</li></ul>
<b>MATERIALS</b>	Polycarbonate   Silicone
<b>CLEANING, USE &amp; STORAGE CONDITIONS</b>	<ul style="list-style-type: none"><li>• Cleaning: Warm water and mild detergent (neutral pH 6 - 8). Do not use solvents (turpentine or acetone), hot water, bleaching or chemical agents</li><li>• Storage conditions: +10°C to +30°C (+50°F to +86°F) at &lt; 75% relative humidity</li><li>• Store away from direct sunlight, grease and oil</li></ul>
<b>APPLICATIONS</b>	<p>CleanSpace HALO Cleaning and Storage Plug is suitable for use during and after cleaning processes, decontamination, and storage.</p> <p>The CleanSpace HALO Cleaning and Storage Plug is designed to be inserted into both bellows and filter air inlet of the CleanSpace HALO / CleanSpace WORK Power Unit.</p>
<b>LIMITATIONS</b>	<p>The plug is not a filter. Do not attempt to run the respirator with the plug fitted. After cleaning, allow Power Unit to air dry prior to inserting cleaning and storage plug for storage.</p>
<b>EDUCATION</b>	Contact <a href="mailto:ordersUSA@cleanspacetechnology.com">ordersUSA@cleanspacetechnology.com</a> or visit the <a href="https://www.cleanspacetechnology.com">CleanSpace website</a> for details.
<b>WARRANTY</b>	<p>Visit <a href="https://www.cleanspacetechnology.com/warranty">www.cleanspacetechnology.com/warranty</a> to read important Warranty information and complete the Warranty Registration Online Form for CleanSpace Products. Consumable parts are warranted up to the point of use, provided they have been stored correctly and are within their expiry date.</p> <p>CleanSpace recommends filters are replaced regularly.</p>

North America (Toll Free): +1 877-256-7723 (1877CLNSPCE)  
[ordersUSA@cleanspacetechnology.com](mailto:ordersUSA@cleanspacetechnology.com)

*Uncontrolled if downloaded or printed from the internet.*

[WWW.CLEANSPACEUSA.COM](https://www.cleanspaceusa.com)

