



# UKCA CERTIFICATE

## for Type Examination (Module B)



<b>Certificate Number:</b>	<b>8509-B-05-22-0061</b>	<b>Amendment:</b>	<b>N/A</b>
<b>Reason for amendment:</b>	N/A		

The following PPE type complies with the applicable Essential Health and Safety Requirements in Annex II of Regulation 2016/425 on personal protective equipment, as amended to apply in GB.

<b>PPE Category:</b>	III
<b>PPE Type:</b>	Power assisted filtering device incorporating full face mask or half masks
<b>Trademark:</b>	CleanSpace ULTRA
<b>Model:</b>	ULTRA
<b>Reference:</b>	CleanSpace ULTRA
<b>Description:</b>	Power assisted filtering device class TM3 P R SL, incorporating a half mask in silicon available in three sizes: S, M and L (references CST1014, CST1015 and CST1016) equipped with one exhalation valve in silicon (housing in polycarbonate), or a full-face mask in silicone available in two sizes : S and M/L (references CST1017 and CST1018) equipped with one exhalation valve in silicon and two inhalation valves in silicon. The battery and the turbo unit reference CST1012 are worn behind the neck. The turbo unit supplies filtered ambient air in demand to the facepiece by a breathing hose with special thread connection adjustable by adjusted clamp. The turbo unit can be equipped either with a standard particles filter reference CST1005 or an extended capacity particles filters reference CST1004 (housing in ABS), both in microfiber glass. The full-face mask is adjustable by a five points head straps in polyethylene and the half mask is adjustable by a three-point head strap in polyethylene.
<b>Issued To:</b>	CleanSpace Technology Pty Ltd - Unit 5 - 39 Herbert Street - ST LEONARDS - NSW 2065 – Australia
<b>Designated Test Standard:</b>	EN 12942:1998 + A1:2002 + A2:2008
<b>Performance Level/Protection Class:</b>	TM3 P R SL
<b>PARC Report Number:</b>	4U00042.2
<b>Date Signed:</b>	04/05/2023
<b>Date of First Issue:</b>	04/05/2023
<b>Date of Renewal:</b>	First edition
<b>Date of Expiry:</b>	04/05/2028

<b>Signatory:</b>	
-------------------	--

### Notes:

1. Issued in accordance with Regulation 2016/425 on personal protective equipment, as amended to apply in GB, and the edition of PARC PPE-R Module B Certification Scheme in force at the issue date of this certificate.
2. For category III PPE, this certificate shall only be used in conjunction with one of the conformity assessment procedures referred to in point (c) of Article 19.
3. The manufacturer shall inform Product Assessment and Reliability Centre Ltd of all modifications to the approved type and all modifications of the Technical Documentation that may affect the conformity of the PPE with the applicable essential health and safety requirements or the conditions for validity of this certificate (Article 7.2, Annex V).

Product Assessment and Reliability Centre Ltd, Unit 4 Alverdiscott Road Industrial Estate, Bideford, Devon, EX39 4LQ, UK

+44 (0) 1237 421255

ukca@parctest.co.uk

[www.parctest.co.uk/](http://www.parctest.co.uk/)