

PROTECTION NEVER LOOKED SO GOOD
WITH THE LATEST MODERN
PAPR TECHNOLOGY



with
AirSensit[®]
TECHNOLOGY

CleanSpace HALO

Approved PAPR with
HEPA filtration; Source control
Lightweight (11lb) and compact
Comfortable, soft silicone mask

Fresh airflow, no fogging and
easy to wear over long periods

**No hoses, belts or waist
mounted battery packs**

Reliable, fast cleaning,
disinfection and sterilization

Reusable and cost effective

Easy to clean, few parts and water
tolerant (IP Rated 66)

Compatible with Steri-Plus
(Source control) and BioHood
System (face/head protection)



LEARN MORE

CleanSpace HALO

Specially designed for clinical care – medical device grade materials, water tolerant and validated for a wide range of cleaning, disinfection and sterilization methods

RESPIRATORS AND MASKS

CS3020	CleanSpace HALO Power System (exc mask)		CS3004	CleanSpace HALO Half Mask MEDIUM
			CS3005	CleanSpace HALO Half Mask LARGE
			CS3006	CleanSpace HALO Full Face Mask SMALL
CS3003	CleanSpace HALO Half Mask SMALL		CS3007	CleanSpace HALO Full Face Mask MEDIUM/LARGE

FILTERS (HEPA), PRE-FILTERS AND SOURCE CONTROL

CS3022	CleanSpace HALO Particulate Filter HEPA (pk 3)		PAF-0058	CleanSpace Particulate Pre-Filter and protective sleeve (pk 20)
CS3039	CleanSpace Steri-Plus Source Control Filters (pk 20) (excl case)		CS3025	CleanSpace HALO BIO HEPA Hood Filter (pk 3)
CS3038	CleanSpace Steri-Plus Source Control Case (excl filters)		CS3027	CleanSpace HALO BIO Hood Exhalation Filter 95 (pk 15)

ACCESSORIES

CS3011	CleanSpace HALO Cleaning & Storage Plug Set		CS3015	CleanSpace HALO Mask ID Badges (Pk 10)
CS3014	CleanSpace Docking Station - Charging & Storage for 8 units		CS3016	CleanSpace Mask Decal (Pk 5)
CS3040	CleanSpace HALO Back Clip for Head Harness		CS3024	CleanSpace BIO Hood (single) compatible with half mask
PAF-0025	CleanSpace Half Mask Adaptor for Portacount Quantitative Fit Test		CS3026	CleanSpace BIO Hood Half Mask Exhalation Adaptor (Pk 5)
PAF-1024	CleanSpace Car Charger		PAF-1017	CleanSpace Full Face Mask Spectacle Kit
CS3013	CleanSpace HALO Unit ID Panels (pk 6)		PAF-1018	CleanSpace Full Face Mask Tear off Visor Anti-Scratch Protectors (pk 10)

CS3020
CleanSpace HALO Power System
(exc mask)

